



Niti-S Enteral Colonic Covered Stent

Head-type



Features

- **Structure : "Head-type" with silicone covering**
 - Fully covered membrane prevents the risk of tumor in-growth
 - Removal suture at the proximal end allows easy and safe removal
- **Radiopaque marker :** 4(four) at the both ends & 2(two) at the middle

Ordering Information

TTS(Through The Scope)

► Full covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806F	18	6	10.5	220
CST1808F		8		
CST1810F		10		
CST1812F		12		
CST1814F		14		
CST1815F		15		
CST2006F	20	6	10.5	220
CST2008F		8		
CST2010F		10		
CST2012F		12		
CST2014F		14		
CST2015F		15		


Coding

CST-----F : S-silicone T-TTS F-full covered

OTW(Over The Wire)

► Full covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806F	18	6	16	70
CS1808F		8		
CS1810F		10		
CS1812F		12		
CS1814F		14		
CS1815F		15		
CS2006F	20	6	20	
CS2008F		8		
CS2010F		10		
CS2012F		12		
CS2014F		14		
CS2015F		15		
CS2206F	22	6	20	
CS2208F		8		
CS2210F		10		
CS2212F		12		
CS2214F		14		
CS2215F		15		
CS2406F	24	6	22	
CS2408F		8		
CS2410F		10		
CS2412F		12		
CS2414F		14		
CS2415F		15		
CS2606F	26	6	22	
CS2608F		8		
CS2610F		10		
CS2612F		12		
CS2614F		14		
CS2615F		15		
CS2806F	28	6	22	
CS2808F		8		
CS2810F		10		
CS2812F		12		
CS2814F		14		
CS2815F		15		



Relative Studies

- "From iatrogenic digestive perforation to complete anastomotic disunion: endoscopic stenting as a new concept of 'stent-guided regeneration and re-epithelialization'"

- Laila Amrani: GASTROINTESTINAL ENDOSCOPY Volume 69, No. 7: 2009

Ordering Information

TTS(Through The Scope)

OTW(Over The Wire)

► Both ends head part bare

► Both ends head part bare



Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806B	18	6	10.5	220
CST1808B		8		
CST1810B		10		
CST1812B		12		
CST1814B		14		
CST1815B		15		
CST2006B	20	6	10.5	220
CST2008B		8		
CST2010B		10		
CST2012B		12		
CST2014B		14		
CST2015B		15		

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806B	18	6	16	70
CS1808B		8		
CS1810B		10		
CS1812B		12		
CS1814B		14		
CS1815B		15		
CS2006B	20	6	20	
CS2008B		8		
CS2010B		10		
CS2012B		12		
CS2014B		14		
CS2015B		15		
CS2206B	22	6	20	
CS2208B		8		
CS2210B		10		
CS2212B		12		
CS2214B		14		
CS2215B		15		
CS2406B	24	6	22	
CS2408B		8		
CS2410B		10		
CS2412B		12		
CS2414B		14		
CS2415B		15		
CS2606B	26	6	22	
CS2608B		8		
CS2610B		10		
CS2612B		12		
CS2614B		14		
CS2615B		15		
CS2806B	28	6	22	
CS2808B		8		
CS2810B		10		
CS2812B		12		
CS2814B		14		
CS2815B		15		

Coding

CST----B : S-silicone T-TTS B-both ends head part bare

Ordering Information



TTS(Through The Scope)

OTW(Over The Wire)

► Distal end head part bare

► Distal end head part bare

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806H	18	6	10.5	220
CST1808H		8		
CST1810H		10		
CST1812H		12		
CST1814H		14		
CST1815H		15		
CST2006H	20	6		
CST2008H		8		
CST2010H		10		
CST2012H		12		
CST2014H		14		
CST2015H		15		


Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806H	18	6	16	70
CS1808H		8		
CS1810H		10		
CS1812H		12		
CS1814H		14		
CS1815H		15		
CS2006H	20	6	20	
CS2008H		8		
CS2010H		10		
CS2012H		12		
CS2014H		14		
CS2015H		15		
CS2206H	22	6	20	
CS2208H		8		
CS2210H		10		
CS2212H		12		
CS2214H		14		
CS2215H		15		
CS2406H	24	6	22	
CS2408H		8		
CS2410H		10		
CS2412H		12		
CS2414H		14		
CS2415H		15		
CS2606H	26	6	22	
CS2608H		8		
CS2610H		10		
CS2612H		12		
CS2614H		14		
CS2615H		15		
CS2806H	28	6	22	
CS2808H		8		
CS2810H		10		
CS2812H		12		
CS2814H		14		
CS2815H		15		

Coding

CST----H : S-silicone T-TTS H-distal end head part bare


Niti-S Enteral Colonic Uncovered Stent

D-type




Features

- **Structure: "D-type"**
 - Excellent conformability facilitates immediate and continuous wall opposition
 - Minimum foreshortening rate for accurate placement of the stent
 - Preoperative bridge to surgery for obstructive colorectal cancer
 - Palliative treatment for obstructing primary left-sided colorectal cancer
- **Radiopaque marker:** 4(four) at the both ends & 2(two) at the middle




Ordering Information



TTS(Through The Scope)					OTW(Over The Wire)						
Code	Stent		Delivery		Code	Stent		Delivery			
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		
CDT1806	18	6	10	220	CD1806	18	6	16	70		
CDT1808		8			CD1808		8				
CDT1810		10			CD1810		10				
CDT1812		12			CD1812		12				
CDT1814		14			CD1814		14				
CDT1815		15			CD1815		15				
CDT2006	20	6			CD2006	20	6	16			
CDT2008		8			CD2008		8				
CDT2010		10			CD2010		10				
CDT2012		12			CD2012		12				
CDT2014		14			CD2014		14				
CDT2015		15			CD2015		15				
CDT2206	22	6			CD2206	22	6	16			
CDT2208		8			CD2208		8				
CDT2210		10			CD2210		10				
CDT2212		12			CD2212		12				
CDT2214		14			CD2214		14				
CDT2215		15			CD2215		15				
CDT2406	24	6			CD2406	24	6	16			
CDT2408		8			CD2408		8				
CDT2410		10			CD2410		10				
CDT2412		12			CD2412		12				
CDT2414		14			CD2414		14				
CDT2415		15			CD2415		15				
Coding CDT---- : D -type structure T -TTS					CD2606	26	6	18			
					CD2608		8				
					CD2610		10				
					CD2612		12				
					CD2614		14				
					CD2615		15				
					CD2806	28	6	18			
					CD2808		8				
					CD2810		10				
					CD2812		12				
					CD2814		14				
					CD2815		15				
					CD3006	30	6	18			
					CD3008		8				
					CD3010		10				
					CD3012		12				
					CD3014		14				
					CD3015		15				


Niti-S Enteral Covered Stent


Comvi-type



Features

- **Structure: "Comvi-type"**
 - Biocompatible PTFE membrane tube is held between inside and outside D-type stent bodies, with each of overlapped layer of the PTFE membrane tube and the inside and outside stent bodies being integrating into a single structure based on D-type
 - Ideal combination of radial force and axial force to maintain full luminal patency in torturous anatomy
- **Radiopaque marker:** 4(four) at the both covered part ends





Ordering Information

TTS(Through The Scope)					OTW(Over The Wire)				
Code	Stent		Delivery		Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CCT1806BA	18	6	10.5	220	CC1806BA	18	6	14	70
CCT1808BA		8			CC1808BA		8		
CCT1810BA		10			CC1810BA		10		
CCT1812BA		12			CC1812BA		12		
CCT2006BA	20	6			CC2006BA	20	6	14	
CCT2008BA		8			CC2008BA		8		
CCT2010BA		10			CC2010BA		10		
CCT2012BA		12			CC2012BA		12		
CCT2206BA	22	6			CC2206BA	22	6	16	
CCT2208BA		8			CC2208BA		8		
CCT2210BA		10			CC2210BA		10		
CCT2212BA		12			CC2212BA		12		
Coding CCT---- BA : C-comvi T-TTS BA-both ends 15mm bare					CC2406BA	24	6	16	
					CC2408BA		8		
					CC2410BA		10		
					CC2412BA		12		
					CC2606BA	26	6	18	
					CC2608BA		8		
					CC2610BA		10		
					CC2612BA		12		
					CC2806BA	28	6	18	
					CC2808BA		8		
					CC2810BA		10		
					CC2812BA		12		

Relative Studies

- "Comparison of efficacies between stents for malignant colorectal obstruction: a randomized, prospective study"

- Semi Park: GASTROINTESTINAL ENDOSCOPY Volume 72, No. 2 : 2010